

Amendments to the Specification:

Please insert Tables 1, 2, and 3 after paragraph [0037] on page 11 of the originally filed Specification as follows:

[0037] System 208 also includes an optional detector 216, which may for example be a CCD camera or other suitable detector for a selected characteristic of the patient's skin. The output from detector 216 is applied to a control 218, which is typically a suitably programmed microprocessor, but may be special purpose hardware or a hybrid of hardware and software. Control 218 controls both the turning on and turning off of source 210 and may also control the power profile of the radiation. Control 218 is also applied to optical system 212 to for example control focus depth for the optical system and to control the portion or portions 214 to which radiation is being focused/concentrated at any given time, for example by controlling scanning by the optical system and/or the beam radiating therefrom. Finally, controls 218 are applied to cooling element 215 to control both the skin temperature above the volume V and the cooling duration, both for precooling and during an irradiation.

Table 1.

Depth of damage , μm	Wavelength range, μm		NA range		Pulse width range, s
	broad	preferred	broad	preferred	
50-200	400 – 1880 & 2050-2350	800-1850 & 2100-2300	<3	0.2 – 1	<2
200-300	500-1880 & 2050-2350	800-1850 & 2150-2300	<3	0.2 – 1	<10
300-500	600-1380 & 1520-1850 & 2150-2260	900-1300 & 1550-1820 & 2150-2250	<2	0.2 – 1	<60
500-1000	600-1370 & 1600-1820	900-1250 & 1650-1750	<2	0.2 – 0.6	<120
1000-2000	670-1350 & 1650-1780	900-1230	<1.5	0.2 – 0.6	<120
2000-5000	800-1300	1050-1220	<1	0.2 – 0.4	<300

Table 2

Depth of damage, μm	Diameter of damage, μm	Wave-length μm	NA	Pulse width, ms	Energy, J	Focusing depth, μm
300	50-100	2.2	0.3-0.5	<10	>0.00015	400-600
300	50-100	1.7	0.3-0.5	<10	>0.0007	400-600
300	50-100	1.3	0.3-0.5	<10	>0.003	400-600
300	50-100	1.54	0.3-0.5	<10	>0.0003	400-600
300	50-100	1.208	0.4-1	<10	>0.016	400-600
300	50-200	0.92	0.4-1	<10	>0.15	400-600
1000	50-200	1.7	0.3-0.4	<100	>0.01	1100-2000
1000	50-200	1.54	0.4	<100	>0.008	1100-2000
1000	50-200	1.3	0.4	<100	>0.1	1100-2000
1000	50-200	1.208	0.4	<100	>0.4	1100-2000

Table 3

Depth of damage, μm	Diameter of damage, μm	Wavelength, μm	NA	Pulse width, ms	Power, W	Focusing depth, μm
500-1000	200-1000	2.2	0.3-0.5	>100	>0.5	600-1500
500-1000	200-1000	1.7	0.3-0.5	>100	>1.5	600-2000
500-1000	200-1000	1.208	0.3-0.6	>3000	>1.0	600-2000
500-1000	400-1200	0.92	0.3-0.6	>3000	>25.0	600-2000
2000-3500	1000-2000	1.208	0.3-0.4	>10000	>1.5	4000-6000